

Search Details

Job Title/ Assay Name : parapsilosis(in)	Amplicon Sequence: GCCTGTTTGAGCGTCATTTCTCCCTCAAACCCTCGGGTTTGGTGTTG AGCGATACGCTGGGTTTGCTTGAAAGAAAGGCGGAGTATAAACTAA TGGATAGGTTTTTTCCACTCATTTGGTACAAACTCCAAAACCTTCTTCC AAATTCGACCTCAAATCAGGTAGGACTACCCGCTGAACTTAAGCAT ATCAATAAGCGGAGG
Effective Date : 2025-11-05	Forward: GCCTGTTTGAGCGTCRTTTC
Database : core_nt	Reverse: CCTCCGCTTATTGATATGCTTA
Alignment algorithm : blastN	Probe: ACTCCGCCTTTCTTTCAAGCAAACCC
Exclusive Organism:	Inclusive Organism: Candida parapsilosis (5480)
RID: GNVJJACU014	Total Result: 1947
App Version: 1.0	

No of Accessions	Accession Details	Oligo	Mutation	BDS Analysis Performance Prediction
Candida parapsilosis/5480 Analyzed Accessions: 1235				
598	598	NA	No Mutations	 Detected
6	KU950722 KY657567 OK303414 MN747310 KT191105 MN747311	Reverse 	C1-	 Detected
34	34	Reverse 	C1-, C2-	 Detected
5	JN989535 EU564204 MZ020644 JN989522 KC777378	Reverse 	T11A	 Detected
8	LC389996 LC389120 LC389985 MT814290 OM523857 OM523858 LC390067 LC389984	Reverse 	T8C	 Detected
2	AB109223 PP733852	Reverse 	T3C	 Detected
3	ON181974 MZ769005 MN747315	Reverse 	C7T	 Detected
3	KM113867 KU052789 PP977577	Reverse 	G6C	 Detected

No of Accessions	Accession Details	Oligo	Mutation	BDS Analysis Performance Prediction
2	MW819639 OM523851	Reverse   	C7G	   Detected
1	KJ706895	Probe   	C21G	   Detected
3	OR287171 MW703694 KJ706736	Reverse   	A10T	   Detected
26	26	Reverse   	C1-, C2-, T3-	   Detected
26	26	Reverse   	+7C	   Detected
27	27	Reverse   	+6G	   Detected
1	EU484061	Reverse   	+6N	   Detected
4	OP895150 KX652405 MK267727 GQ280288	Reverse   	+6T	   Detected
6	KP131782 KP131764 KP674489 JN989536 KJ734198 MW703702	Reverse   	+9T	   Detected
2	AB109224 AB109231	Reverse   	+5N	   Detected
3	JN989524 JN989516 JN989532	Reverse   	+10A	   Detected
4	LC390139 LC389118 KP674743 LC390144	Reverse   	+6A	   Detected
16	16	Reverse   	C1-, C2-, T3-, C4-	   Detected
1	FR774537	Forward   	G1-, C2-, C3-, T4-	   Detected
5	LC390062 LC389954 LC389972 KM113926 OM523853	Reverse   	T9-	   Detected
3	KM113884 KM113880 KM113895	Reverse   	G6-	   Detected
3	KT067766 OQ184042 PV130509	Reverse   	+6G C	   Detected
1	MG599213	Reverse   	C1-, +9T	   Detected

No of Accessions	Accession Details	Oligo	Mutation	BDS Analysis Performance Prediction
1	GQ280301	Reverse   	+5C G	   Detected
1	OR091288	Reverse   	T12-	   Detected
1	KP675355	Reverse   	A10-	   Detected
1	MT729884	Forward  	T14C	  Detected
14	EU871509 LC389122 OW985494 OW986219 LC389810 KP783504 LC389139 OW983289 OK303405 OK317693 MW709969 LC389787 OW984120 OL687495	Reverse   	C1-, C2-, T3-, C4-, C5-	   Detected
10	OL774548 MH299862 MK392005 MW819637 KP674744 KM113857 KM113990 KM113890 AB109235 OM523863	Reverse   	C7G, T8C	   Detected
1	OM959548	Reverse   	A10T, T11A	   Detected
1	LC390097	Reverse   	C7A, T8C	   Detected
3	PV363933 OM237275 PV363932	Reverse   	T8C, A10T	   Detected
2	LC390068 KM113862	Reverse   	C7T, T8C	   Detected
4	MT815460 JQ993383 LC389986 AB109288	Reverse   	C1-, C2-, T3-, T8C	   Detected
1	MN450875	Reverse   	C1-, T9-	   Detected
3	OM959620 OM959541 OM959655	Reverse   	C1-, C2-, T3-, A10T	   Detected
1	KP131759	Reverse   	C1-, C2-, T3-, G6C	   Detected
1	MN450869	Forward   	T8-	   Detected
		Reverse   	C1-	

No of Accessions	Accession Details	Oligo	Mutation	BDS Analysis Performance Prediction
1	MG599241	Reverse   	C1-, C2-, T3-, T11A	   Detected
1	KY075668	Reverse   	C1-, C2-, +7C	   Detected
1	PQ248795	Reverse   	+6N, T8C	   Detected
19	19	Reverse   	C1-, C2-, T3-, C4-, C5-, G6-	   Detected
1	KP131758	Reverse   	A10T, +12T	   Detected
1	OM523859	Reverse   	C7G, +7G	   Detected
1	OM523862	Reverse   	C7G, +7A	   Detected
1	LC389119	Reverse   	G6A, T9-	   Detected
40	40	Reverse   	C1-, C2-, T3-, C4-, C5-, G6-, C7-	   Detected
2	AB109286 AB109282	Reverse   	T3C, T9-	   Detected
1	PV363916	Reverse   	C1-, C2-, T8C, A10T	   Detected
1	AB109292	Reverse   	T3C, C7G, T8C	   Detected
8	AB109233 AB109290 AB109229 AB109291 AB109230 AB109234 KM113868 AB109281	Reverse   	C1-, C2-, T3-, T9-	   Detected
1	KP674827	Reverse   	G6C, C7G, A10T	   Detected
5	PQ206277 OP618190 PQ206278 KP674736 KP674519	Reverse   	C1-, C2-, T3-, A10-	   Detected
1	PP844419	Probe    Reverse   	G25A C7A	   Detected
3	KP674908 MH166866 KP674969	Reverse   	C7G, T8C, A10T	   Detected
1	KP674652	Reverse   	C7A, T8C, A10T	   Detected
1	EF195357	Reverse 	C1-, C2-, T3-, T12-	 Detected
44	44	Reverse   	C1-, C2-, T3-, C4-, C5-, G6-, C7-, T8-	   Detected

No of Accessions	Accession Details	Oligo	Mutation	BDS Analysis Performance Prediction
3	MF462176 MK503153 KP675363	Reverse   	C1-, C2-, T3-, T8C, A10T	   Detected
5	OP618179 OP618188 OP618187 OP618193 MT153735	Reverse   	+6G, T12-	   Detected
1	PV074065	Reverse   	C1-, C2-, T3-, C7G, A10T	   Detected
1	OM681576	Reverse   	+5S, T8C, A10T	   Detected
1	OM681577	Reverse   	C1-, C2-, T3-, C4-, C5-, G6-, A10T	   Detected
3	MW709971 KP675635 KP675616	Reverse   	+6G, T8C, A10T	   Detected
1	OR784796	Probe   	G20A	   Detected
24	24	Reverse   	C1-, C2-, T3-, C4-, C5-, G6-, C7-, T8-, T9-	   Detected
10	EF193066 KP674607 MT729857 MF380666 EF198017 LC790897 MT729856 OP618194 MK392002 MT729855	Reverse   	+6G G, A10-	   Detected
2	KP674651 KP675285	Reverse   	+5C G, A10-	   Detected
1	OP618189	Reverse   	+5C G, T12-	   Detected
1	AB109284	Reverse   	T3C, A10-, T11-	   Detected
8	JN989527 JN989519 JN989525 JN989517 JN989537 KX652403 MH445556 JN989533	Reverse   	C1-, C2-, T3-, C4-, C5-, G6-, C7-, T11A	   Detected
1	PP033709	Reverse   	T3C, T8C, A10-	   Detected
1	PV363915	Reverse   	C1-, C2-, T3-, C4-, C5-, G6-, C7-, A10T	   Detected
1	QQ876043	Reverse 	C1-, C2-, T3-, C4-, C5-, T12-	 Detected
6	KC422431 OR507588 KP675663	Reverse   	C1-, C2-, T3-, C4-, C5-, A10-	   Detected

No of Accessions	Accession Details	Oligo	Mutation	BDS Analysis Performance Prediction
	KP204942 KU961982 KP674958			
1	EF191046	Reverse   	T3C, G6C, A10-	   Detected
1	KP131779	Reverse   	+6G, A10T, +11A	   Detected
71	71	Reverse   	C1-, C2-, T3-, C4-, C5-, G6-, C7-, T8-, T9-, A10-	   Detected
1	MT729838	Reverse 	C1-, C2-, T3-, G6C, T12-	 Detected
1	KJ706418	Probe   	A24G	   Detected
		Reverse   	T3C	
5	OM959604 PQ796639 EF187252 OM959551 PV363934	Reverse 	C1-, C2-, T3-, C4-, C5-, G6-, T12-	 Detected
1	OP603449	Reverse   	C1-, C2-, T3-, G6C, A10T, T11A	   Detected
60	60	Reverse   	C1-, C2-, T3-, C4-, C5-, G6-, C7-, T8-, T9-, A10-, T11-	   Detected
1	OR820540	Probe   	T6-	   Detected
		Reverse   	C1-, C2-	
1	MF767673	Probe   	G1A, C8T	   Detected
1	OM523860	Probe   	G22C	   Detected
		Reverse   	G6T	
1	KY560334	Probe  	G15A	   Detected
		Reverse   	C1-, C2-, T3-, C4-, C5-, G6-, C7-, T8-, T9-, A10-	
1	PP195927	Probe  	+14AG	   Detected
		Reverse   	C1-, C2-, T3-, C4-, C5-, G6-, C7-, T8-, T9-, A10-, T11-	
1	EU258672	Forward   	T8-	   Detected
		Probe  	T9C	
1	JX463241	Probe   	+3G	   Detected
1	KU745529	Forward   	G11-	   Detected
		Reverse   	C1-, C2-, T3-	
1	KY075671	Probe   	T6A	   Detected
		Reverse   	C1-, C2-, T3-, C4-, C5-, G6-, C7-, T8-	
1	JN989526	Reverse   	T11A, A22T	   Not Detected
1	KU961981	Reverse   	C1-, C2-, T3-, C4-, C5-, C19G	   Not Detected

No of Accessions	Accession Details	Oligo	Mutation	BDS Analysis Performance Prediction
1	JX463245	Reverse XXX	T20C	XXX Not Detected
1	EF192217	Reverse X	C1-, T8C, A10C, +12T	X Not Detected
1	KJ451700	Reverse X	C1-, C2-, T9C, T11C, T12C	X Not Detected
1	MW132094	Forward XXX	C3G, A16G	XXX Not Detected
		Reverse XXX	T20C	
1	MF510442	Reverse XXX	G13C, T15G, T17G, T21G	XXX Not Detected
2	MT876148 MT876137	Reverse X	C1-, C2-, T3-, C4-, C5-, G6-, C7-, T8-, T9-, A10-, G13T	X Not Detected
1	MT729885	Reverse XXX	C4G, C7G, T11C, T17A, C19G	XXX Not Detected
1	OK090870	Reverse X	C1-, C2-, T3-, G6C, A10T, A14C	X Not Detected
1	KY560331	Forward XXX	A16T	X Not Detected
		Probe ✓✓	G15A	
		Reverse ✓✓✓	C1-, C2-, T3-, C4-, C5-, G6-, C7-, T8-	
1	EU564209	Forward XXX	C20-	XXX Not Detected
1	OR243746	Reverse XXX	C1-, C2-, T3-, C4-, C5-, G6-, C7-, C19G, T20G	XXX Not Detected
1	KC182136	Probe ✓✓✓	C21G	XX Not Detected
		Reverse XXX	C1-, +6G, C19G, T21A	
1	MW703799	Reverse XXX	T12G, A16C, T17C, T20G	XXX Not Detected
2	MT322428 MT322429	Forward ✓✓✓	T8C	X Not Detected
		Probe ✓✓	G11C	
		Reverse XXX	A10G, T11A, T21G	
1	OQ520002	Reverse XXX	T9-, C19T	XXX Not Detected
1	MF767641	Forward XXX	A16G	XXX Not Detected
7	PP800201 PP800203 MT322430 PP800202 PP800204 PP800206 PP800205	Forward ✓✓✓	T8C	X Not Detected
		Probe ✓✓	G11C	
		Reverse XXX	C1-, C2-, T8G, A10G, T11A, T21G	
1	JX463259	Probe XXX	+3G, +14A, +23G	XXX Not Detected
1	JX406279	Forward XX	G9C, C12A, C15G	XX Not Detected
		Probe ✓✓✓	G3A	
		Reverse ✓✓✓	T11A	

No of Accessions	Accession Details	Oligo	Mutation	BDS Analysis Performance Prediction
1	MT729948	Forward XXX	G11C, T14C, A16G	XXX Not Detected
		Probe XXX	A14G, C21G, G25C, T26G	
		Reverse X	C1-, C2-, T3-, C4-, C5-, G6-, C7-, T8-, T9-, A10-, T11-	
1	MK886556	Forward XXX	T19G	XX Not Detected
		Probe XXX	T10C, A14-, T26-, +22A	
		Reverse ✓✓✓	C7T	
3	MT605316 MT605317 MT605315	Forward ✓✓✓	T8C	XX Not Detected
		Probe XX	G11C, T26C	
		Reverse XXX	C1-, C2-, T8G, A10G, T11A, T21G	
1	OK090872	Probe ✓✓✓	G1T	XX Not Detected
		Reverse XXX	C1-, C2-, T3-, C4-, C5-, G6-, C7-, T8-, T9-, A10-, G13A, A16G, C19G, A22G	
1	MW506001	Reverse XXX	C1-, C5G, +15G, G18A, T20C, T21C	XXX Not Detected
1	JN942629	Probe XX	+15A, +25A	XX Not Detected
		Reverse X	C1-, C2-, T3-, C4-, C5-, G6-, C7-, T8-, T9-, A10-	
1	KY106657	Forward X	G1-, C2-, C3-, T4-, G5-, T6-, T7-, T8-, G9-, A10-, G11-	X limited coverage (excluded)
59	59	Reverse X	C1-, C2-, T3-, C4-, C5-, G6-, C7-, T8-, T9-, A10-, T11-, T12-	X limited coverage (excluded)
12	OL334739 PV739597 PP858936 PP858935 MH782053 PV384124 PQ570470 OR199906 PP218380 MT378269 PP204162 MW940754	Reverse X	C1-, C2-, T3-, C4-, C5-, G6-, C7-, T8-, T9-, A10-, T11-, T12-, G13-	X limited coverage (excluded)
23	23	Reverse X	C1-, C2-, T3-, C4-, C5-, G6-, C7-, T8-, T9-, A10-, T11-, T12-, G13-, A14-	X limited coverage (excluded)
19	19	Reverse XX	C1-, C2-, T3-, C4-, C5-, G6-, C7-, T8-, T9-, A10-, T11-, T12-, G13-, A14-, T15-	XX limited coverage (excluded)
1	KC556814	Reverse XXX	C1-, C2-, T3-, C4-, C5-, G6-, C7-, T8-, T9-, A10-, T11-, T12-, G13-, A22G	XXX limited coverage (excluded)
6	PV416740 PV053617 KP131791	Reverse XXX	C1-, C2-, T3-, C4-, C5-, G6-, C7-, T8-, T9-, A10-, T11-, T12-, G13-, A14-, T15-, A16-	XXX limited coverage (excluded)

No of Accessions	Accession Details	Oligo	Mutation	BDS Analysis Performance Prediction
	KT191118 PP229797 KU987850			
429	429	Reverse XXX	Reverse Missing	XXX limited coverage (excluded)
1	KY560339	Probe ✓✓	G15A	✗ limited coverage (excluded)
		Reverse ✗	C1-, C2-, T3-, C4-, C5-, G6-, C7-, T8-, T9-, A10-, T11-, T12-	
1	PP195928	Probe ✓✓	+15A	✗ limited coverage (excluded)
		Reverse ✗	C1-, C2-, T3-, C4-, C5-, G6-, C7-, T8-, T9-, A10-, T11-, T12-	
1	KY560336	Probe ✓✓	G15A	✗✗ limited coverage (excluded)
		Reverse XXX	C1-, C2-, T3-, C4-, C5-, G6-, C7-, T8-, T9-, A10-, T11-, T12-, G13-, A14-, T15-, A16-	
1	KY560337	Probe ✓✓	G15A	✗ limited coverage (excluded)
		Reverse ✗✗	C1-, C2-, T3-, C4-, C5-, G6-, C7-, T8-, T9-, A10-, T11-, T12-, G13-, A14-, T15-	
1	KC858280	Probe ✓✓	G11C	✗ limited coverage (excluded)
		Reverse ✗	C1-, C2-, T3-, C4-, C5-, G6-, C7-, T8-, T9-, A10-, T11-, T12-, G13-, A14-	
1	KY560332	Probe ✓✓	G15A	✗✗ limited coverage (excluded)
		Reverse XXX	Reverse Missing	
6	LN890938 MH844406 LT549053 LT549062 KY106663 PV826898	Forward XXX	Forward Missing	XXX limited coverage (excluded)
1	PQ345548	Probe ✓✓	A16-	✗✗ limited coverage (excluded)
		Reverse XXX	Reverse Missing	
2	KY106655 KY106649	Forward XXX	Forward Missing	✗ limited coverage (excluded)
		Reverse ✓✓✓	+6G	
1	MN365727	Forward ✓✓✓	G9T	✗ limited coverage (excluded)
		Reverse XXX	Reverse Missing	
1	KY106658	Forward XXX	Forward Missing	✗✗ limited coverage (excluded)
		Probe ✓✓✓	+21C	
1	PQ363652	Probe ✓✓	A13G	✗✗ limited coverage (excluded)
		Reverse XXX	Reverse Missing	
1	KX129764	Probe XXX	A14-, A18-, G19-	XXX limited coverage (excluded)
		Reverse XXX	Reverse Missing	
1	KY106651	Forward XXX	Forward Missing	✗✗ limited coverage (excluded)
		Probe ✓✓✓	+23G	

No of Accessions	Accession Details	Oligo	Mutation	BDS Analysis Performance Prediction
1	MG496024	Probe   	+20GC	  limited coverage (excluded)
		Reverse   	Reverse Missing	
2	PP838313 PP838316	Probe  	G11-	  limited coverage (excluded)
		Reverse   	Reverse Missing	
1	KP794152	Probe   	G20-	  limited coverage (excluded)
		Reverse   	Reverse Missing	
1	GQ497158	Probe   	G25-	  limited coverage (excluded)
		Reverse   	Reverse Missing	
1	PP838307	Probe  	T10-	  limited coverage (excluded)
		Reverse   	Reverse Missing	
1	PP838311	Probe   	C8-	  limited coverage (excluded)
		Reverse   	Reverse Missing	
1	MT814294	Forward  	T6G, T14G	  Further Analysis
		Probe   	T26G	
		Reverse	T8C, T11A, +14A	
1	MZ566838	Forward   	G11C	 Not Detected
		Probe   	G1T, C21G	
		Reverse   	C1-, C2-, A16G, T17G, G18T, C19A, T20G	
1	LC643454	Probe   	G1A, A16T, A17T, C21A	   Not Detected
1	OM523864	Forward   	C3G, G9C, G11C, G13A, A16G	   Not Detected
		Probe   	G1T, G15A, A16C, A17T, A18C, T26A	
		Reverse   	C7A, T8G, T9A, T20C	
2	MN826506 MW080686	Forward   	A16G	   Not Detected
		Probe   	A24G, G25A, +16CGGTAGTGGT	
1	PV536168	Forward   	A16G	  Not Detected
		Probe   	A24G, G25A, +16CGGTAGTGGT	
		Reverse   	+7C	
1	LC317525	Probe   	G22C, G25C	  limited coverage (excluded)
		Reverse   	Reverse Missing	
1	JF715190	Probe   	G15A, C21G, T26G	   limited coverage (excluded)
		Reverse   	Reverse Missing	
1	MF510440	Probe  	A13C	  limited coverage (excluded)
		Reverse   	Reverse Missing	
1	KU987857	Probe   	G20C	  limited coverage (excluded)
		Reverse   	Reverse Missing	

No of Accessions	Accession Details	Oligo	Mutation	BDS Analysis Performance Prediction
<u>1</u>	KY106675	Forward XXX	Forward Missing	XXX limited coverage (excluded)
		Probe XXX	G1-, G2-, G3-, T4-, T5-	
<u>2</u>	KY106662 KY106664	Forward XXX	Forward Missing	XX limited coverage (excluded)
		Probe XXX	G1-, G2-, G3-, T4-, T5-, T6-, G7-, C8-	
		Reverse ✓✓✓	+6G	
<u>2</u>	KY106660 KY106665	Forward XXX	Forward Missing	XXX limited coverage (excluded)
		Probe XXX	G1-, G2-, G3-, T4-, T5-, T6-, G7-, C8-, T9-, T10-, G11-, A12-, A13-	
<u>1</u>	PQ561357	Forward XXX	Forward Missing	XXX limited coverage (excluded)
		Probe XXX	G1-, G2-, G3-, T4-, T5-, T6-, G7-, C8-, T9-, T10-, G11-, A12-, A13-, A14-, G15-, A16-	
<u>1</u>	KY106666	Forward XXX	Forward Missing	XX limited coverage (excluded)
		Probe XXX	G1-, G2-, G3-, T4-, T5-, T6-, G7-, C8-, T9-, T10-, G11-, A12-, A13-, A14-	
		Reverse ✓✓✓	+6G	
<u>1</u>	AY207079	Forward XXX	G1-, C2-, C3-, T4-, G5-, T6-, T7-, T8-, G9-, A10-, G11-, C12-, G13-, T14-, C15-, A16-, T17-, T18-, T19-	XXX limited coverage (excluded)
		Reverse XXX	Reverse Missing	
<u>1</u>	MK886565	Forward ✓✓✓	T8C	X limited coverage (excluded)
		Reverse XXX	Reverse Missing	
<u>1</u>	HF934031	Forward XXX	Forward Missing	XXX limited coverage (excluded)
		Probe XXX	G1-, G2-, G3-, T4-, T5-, T6-, G7-, C8-, T9-, T10-, G11-, A12-, A13-, A14-, G15-, A16-, A17-, A18-, G19-	
<u>1</u>	PQ536622	Forward XXX	Forward Missing	XXX limited coverage (excluded)
		Probe XXX	G1-, G2-, G3-, T4-, T5-, T6-, G7-, C8-, T9-, T10-, G11-, A12-, A13-, A14-, G15-, A16-, A17-	
<u>1</u>	KY611781	Forward XXX	+16A	XX limited coverage (excluded)
		Probe ✓✓	+15A	
		Reverse XXX	Reverse Missing	
<u>1</u>	KM115147	Forward XXX	G9A, T14G, C15G	XXX limited coverage (excluded)
		Probe XXX	G2T, G20C, G23C	
		Reverse XXX	Reverse Missing	
<u>1</u>	MH854532	Probe ✓✓	+13AG	XX limited coverage (excluded)
		Reverse XXX	Reverse Missing	
<u>20</u>	<u>20</u>	Forward XXX	Forward Missing	XXX limited coverage (excluded)
		Probe XXX	Probe Missing	
<u>1</u>	LC632671	Probe XXX	G19A, C21G, T26G	XXX limited coverage (excluded)
		Reverse XXX	Reverse Missing	
<u>1</u>	MN749538	Forward ✓✓✓	G1-	X limited coverage (excluded)
		Reverse XXX	Reverse Missing	

No of Accessions	Accession Details	Oligo	Mutation	BDS Analysis Performance Prediction
1	MW703698	Probe  	A12G	  limited coverage (excluded)
		Reverse   	Reverse Missing	
1	LC317526	Forward   	Forward Missing	  limited coverage (excluded)
		Probe   	Probe Missing	
		Reverse   	C1-, C2-, T3-, C4-, C5-, G6-, C7-	
1	KY604746	Probe   	G25-, T26-	  limited coverage (excluded)
		Reverse   	Reverse Missing	
2	KM361840 KM361843	Probe   	G22-, G23-, A24-, G25-, T26-	   limited coverage (excluded)
		Reverse   	Reverse Missing	
1	PV687367	Probe   	C21-, G22-, G23-, A24-, G25-, T26-	   limited coverage (excluded)
		Reverse   	Reverse Missing	
1	PP858937	Probe   	G20-, C21-, G22-, G23-, A24-, G25-, T26-	   limited coverage (excluded)
		Reverse   	Reverse Missing	
1	PP504053	Probe   	A18-, G19-, G20-, C21-, G22-, G23-, A24-, G25-, T26-	   limited coverage (excluded)
		Reverse   	Reverse Missing	
1	MT152379	Probe   	G15-, A16-, A17-, A18-, G19-, G20-, C21-, G22-, G23-, A24-, G25-, T26-	   limited coverage (excluded)
		Reverse   	Reverse Missing	
1	MN450870	Forward   	T7C, C12A, T18C	  limited coverage (excluded)
		Probe   	T6C, T26C	
		Reverse   	Reverse Missing	
2	DQ141231 DQ141230	Probe   	G1A, A16T, A17T, C21A	   limited coverage (excluded)
		Reverse   	Reverse Missing	
2	KY611789 MN080443	Probe   	T6-, G7-, C8-, T9-, T10-, G11-, A12-, A13-, A14-, G15-, A16-, A17-, A18-, G19-, G20-, C21-, G22-, G23-, A24-, G25-, T26-	   limited coverage (excluded)
		Reverse   	Reverse Missing	
1	OQ099751	Probe   	G2-, G3-, T4-, T5-, T6-, G7-, C8-, T9-, T10-, G11-, A12-, A13-, A14-, G15-, A16-, A17-, A18-, G19-, G20-, C21-, G22-, G23-, A24-, G25-, T26-	   limited coverage (excluded)
		Reverse   	Reverse Missing	
1	MN309729	Forward  	G11A, C12A	  limited coverage (excluded)
		Probe   	C8G, A14T, A17C, A18T, C21A, G22T, T26C	
		Reverse   	Reverse Missing	
16	16	Probe   	Probe Missing	   limited coverage (excluded)
		Reverse   	Reverse Missing	
3	MZ882284 OK090877 LC713269	Forward   	A16G	   limited coverage (excluded)
		Probe   	A17-, A18-, G19-, G20-, C21-, G22-, G23-, A24-, G25-, T26-	
		Reverse   	Reverse Missing	

No of Accessions	Accession Details	Oligo	Mutation	BDS Analysis Performance Prediction
1	LC534405	Forward 	C15T	 limited coverage (excluded)
		Probe 	Probe Missing	
		Reverse 	Reverse Missing	
1	PV240541	Forward 	Forward Missing	 limited coverage (excluded)
		Probe 	Probe Missing	
		Reverse 	C1-, C2-, T3-, C4-, C5-, G6-, C7-, T8-	
1	KT876508	Forward 	Forward Missing	 limited coverage (excluded)
		Probe 	Probe Missing	
		Reverse 	C1-, C2-, T3-, C4-	
1	JF916558	Probe 	A13-, A14-, G15-, A16-, A17-, A18-, G19-, G20-, C21-, G22-, G23-, A24-, G25-, T26-	 limited coverage (excluded)
		Reverse 	Reverse Missing	
1	OK605288	Forward 	Forward Missing	 limited coverage (excluded)
		Probe 	Probe Missing	
		Reverse 	+4T, T15C	
1	MN365726	Forward 	G1-, C2-, G11C	 limited coverage (excluded)
		Probe 	G1C, T4G, A12-, A13-, A14-, G15-, A16-, A17-, A18-, G19-, G20-, C21-, G22-, G23-, A24-, G25-, T26-	
		Reverse 	Reverse Missing	
1	MZ203560	Forward 	G9A	 limited coverage (excluded)
		Probe 	Probe Missing	
		Reverse 	Reverse Missing	
1	KX268323	Probe 	G1C, G3T, G7-, C8-, T9-, T10-, G11-, A12-, A13-, A14-, G15-, A16-, A17-, A18-, G19-, G20-, C21-, G22-, G23-, A24-, G25-, T26-	 limited coverage (excluded)
		Reverse 	Reverse Missing	
48	48	Forward 	Forward Missing	 limited coverage (excluded)
		Probe 	Probe Missing	
		Reverse 	Reverse Missing	
1	LC617352	Forward 	C15-	 limited coverage (excluded)
		Probe 	Probe Missing	
		Reverse 	Reverse Missing	
2	OR507647 XR_005013741	Forward 	C20-	 limited coverage (excluded)
		Probe 	Probe Missing	
		Reverse 	Reverse Missing	